

Laura Wallace
Sensible Technologies, Inc.
781-939-7437
media@Sensible.com

Mary Kae Marinac
MKM Corporate Communications
978-685-3136
mkm@mkmarinac.com

FOR IMMEDIATE RELEASE

**Sensible Dental Launches New European Version of its Intellifit™
Dental Lab Solutions at IDS**

***Adds German Interface, More Productivity and Production Choice from “Intelligently Open”
Configurations of Industry’s #1 Selling TE CAD/CAM for Partials and More***

WILMINGTON MA March 16, 2011 – [Sensible Dental](#), a division of Sensible, announced that it is launching the newest version of its [Intellifit™](#) TE (Touch-Enabled) Digital Restoration System for European markets at the [34th International Dental Show \(IDS\)](#) in Cologne, Germany, March 22-26. With a new German language interface and the industry’s first virtual touch-enabled articulator so technicians can “feel the fit” as they check occlusion, the new modular configurations of Intellifit are designed expressly for the needs of the European dental lab market. The first and #1 selling system for partials and more, Intellifit – formerly named the Sensible Dental Lab System – now adds support for veneers, enhanced design features and fabrication choices for crown and bridge, as well as increased 3D printing capacity.

Intellifit’s “intelligently open™” platform offers European labs the best of all worlds. The system’s open, modular configurations are ideal for small labs that don’t yet have a CAD/CAM system, as they can choose to purchase just the fixed restoration module, just the removables module, or both. It’s also a great choice for labs that may already have a crown and bridge system and now want to add the industry’s leading solution for digitally designed partials. With Intellifit, European labs can choose to design and fabricate in-house, or outsource to an [Intellifit Production Center](#) for milled, cast or printed parts. Whether labs are large or small, Sensible Dental’s Intellifit solutions help them achieve greater productivity, profitability, and precision to gain a competitive advantage.

“Sensible’s family of 3D design products have been helping the dental industry solve the complex challenges of organic design for the past 10 years,” stated Joerg Bressemer, master dental technician who owns the Bressemer Zahntechnik dental lab. “From modeling many of the world’s most renowned tooth libraries, to helping dental surgeons design custom implants to the new Intellifit Dental Restoration System, technicians have come to respect the range of capabilities the Sensible 3D design system provides. It delivers ultimate flexibility for the technician, from simple cases to the most complex.”

“The Intellifit system is the only digital tool that totally reflects the lab’s needs – with no limitations,” Bressemer added.

Highlights of the new versions of Intellifit on display at IDS include:

- **Enhanced Crown & Bridge Functionality.** The new Intellifit adds the Heraeus Premium® digital tooth library – renowned for its distinctively natural appearance and perfect fit – allowing dental labs even more choice in creating aesthetically pleasing restorations for their doctors. The new Intellifit’s streamlined margin detection also delivers precise margin lines for even the most challenging cases. Multi-die scanning provides accurate 3D scans in less than 20 seconds per unit. And its new touch-enabled, virtual articulator allows technicians to test occlusion of restorations — before they are produced.
- **Support for prepped veneers.** Intellifit eliminates tedious, time-consuming hand-waxing, ensures symmetry and provides consistency, regardless of which technician is doing the work.
- **More Pressing Choices.** Dental labs using Intellifit can design and fabricate full anatomy crowns, bridges and veneers and then press in popular all-ceramic materials from VITA and Ivoclar Vivadent such as Ivoclar Vivadent’s IPS e.max® and IPS Empress® Esthetic.
- **More Milling Choices,** including automatic tool compensation and output for industry standard milling machines. Labs can mill full contour overstructures and substructures in a variety of materials including zirconia and all-ceramic materials through Sensable’s network of authorized milling centers, including Core 3D’s worldwide network
- **Increased Printer Throughput,** with an up to 60% larger printing area, now available on the integrated, high-resolution 3D Systems ProJet™ DP 3000 resin printer, for increased productivity. Labs can print up to 200 parts in a single build. High-resolution resin parts can be used for casting, flasking and pressing and are ideal for fixed and removable restorations.

Sensable will also be demonstrating its Freeform™ 3D Organic Design System, used to create tooth libraries, dental aligners, surgical guides and specialized appliances, and the leading organic modeling system for the design of custom medical implants.

Milling and fabrication centers should visit Sensable to discuss opportunities to become an Intellifit Production Center partner, as a complement to the company’s worldwide network of design users. Sensable is also seeking additional distributors for the Intellifit Digital Restoration System.

Sensable is showcasing these solutions at IDS, Cologne, in Hall 4.1, Stand F-031, beginning Tuesday March 22 through Saturday March 26. For more information visit www.sensabledental.com or contact Sensable at 781.937.8315.

About Sensable

Founded in 1993, Sensable remains the leading developer of touch-enabled solutions and technology that allow users to not only see and hear an on-screen computer application, but to actually “feel” it. With 44 patents granted and over 9,000 systems installed worldwide, Sensable helps people innovate with human touch solutions. The company’s Sensable Dental division provides the Intellifit Digital Restoration System, while Sensable also sells a suite of 3D design solutions that includes its flagship product, Freeform®; and the Phantom® and Omni™ lines of haptic devices, used in surgical simulation and planning, stroke rehabilitation and medical training, as well as a range of research and robotic applications. With an unparalleled commitment to partnering with customers, Sensable brings a human

touch to innovating and implementing customer-centric solutions. Sensable products are available through direct and reseller channels worldwide. www.SensableDental.com.

#

Freeform, Intellifit, Omni, Phantom, "Intelligently Open," Sensable and Sensable Technologies, Inc. are trademarks or registered trademarks of Sensable Technologies, Inc. Other brand and product names are trademarks of their respective holders.